

CE
0086

V7 3/18

Light Bond™ Sealants

**With or Without FLUORIDE
LIGHT CURE SEALANT RESIN**

INDICATIONS FOR USE: Light Bond™ sealants are intended for use as a light cure sealant resin for bonding under an orthodontic adhesive.

INSTRUCTIONS FOR USE:

Dispense one to two drops of Light Bond sealant resin onto a mixing pad. With a brush, apply a thin, uniform layer of sealant to the etched enamel. It is not necessary to cure the sealant at this time; however, a cured sealant can prevent bracket floatation.



PRECAUTIONS: Unpolymerized resins may cause skin sensation in susceptible persons. In case of contact with the skin, wash thoroughly with soap and water.



Flammable



TEMPERATURE LIMITATION: Store at room temperature (15-30 °C)

TO REORDER: Ref# LBS3F – 3cc Light Bond Sealant w/Fluoride, LBS7F – 7cc Light Bond Sealant w/Fluoride, LBS3 – Non-Fluoride, 3cc Light Bond Sealant, LBS7 – Non-Fluoride 7cc Light Bond Sealant, LBFS7F – Light Bond Filled Sealant 7cc w/Fluoride, LBFSF – Light Bond Filled Sealant 3cc w/ Fluoride

WARRANTY: Reliance Orthodontic Products, Inc. recognizes its responsibility to replace products if proven to be defective. Reliance Orthodontic Products, Inc. does not accept liability for any damages or loss, direct or consequential, stemming from the use of or inability to use the products as described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

Rx Only: Federal law restricts this device to sale by or on the order of a dental professional.



Reliance Orthodontic Products, Inc.



1540 West Thorndale Ave., Itasca, IL 60143 U.S.A.
1-800-323-4348 • www.RelianceOrthodontics.com



Scan for more information



EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Australian Sponsor:
Emergo Australia
Level 20, Tower II, Darling Park
201 Sussex St.
Sydney, NSW 2000 Australia