

WaterPik® WP-660 Aquarius Professional Water Flosser

Featuring a new compact and contemporary design, the Aquarius® includes on/off water control on the handle, an LED information panel, and 7 water flosser tips. New pulse-modulation technology provides maximum plaque removal (Floss Mode) and enhanced gum stimulation and improved circulation (Hydro-Pulse Massage Mode). The convenient one-minute timer with a 30-second pacer ensures thorough water flossing of all areas. It includes 3 - Classic Jet Tip, 1 - Plaque Seeker Tip, 1 - Orthodontic Tip, 1 - Pik Pocket Tip, 1 Tooth Brush Tip



Product Features :

- 5 To 90 Psi
- Advanced Pressure Control System With 10 Settings
- Reservoir With Built-In Tip Storage
- Pause Button On Handle
- Reservoir Capacity Of 90 Seconds+Smaller And Quieter Than Previous Models.

WaterPik® Sensonic® Professional Plus Toothbrush SR-3000

Rechargeable
no batteries
needed

Introducing the new **WaterPik® Sensonic® Professional Plus** toothbrush. With state-of-the-art sonic technology..

The advanced brush head design, with extra soft, end-rounded bristles gently targets those hard to reach areas between teeth. It also features a 2 minute brushing timer with 30 second intervals to ensure adequate brushing time throughout the mouth, hygienic brush head covers and a premium travel case.

Naturally whitens teeth and promotes healthy gums

Product Features

- New charging base and deluxe recharge indicator gauge
- Three brush head options to meet individual needs
- Dual speed control on handle: high for superior plaque removal and low for gentle cleaning and gum massage
- Premium Travel Case
- Hygienic Brush Head Covers • 30 second quadrant pacer



Includes • 1 - Standard Brush Head • 1 - Compact Brush Head
• 1 - Interdental Brush Head

WaterPik® FLW-220 Power Flossers

- Clinically proven as effective as string floss and easier to use.
- Cleans and whitens between teeth where plaque forms and stains are toughest to reach.
- 10,000 gentle strokes per minute remove plaque above and below the gum line.
- Flexible nylon tip slips easily between teeth and is gentle enough on gums to use every day.
- Ergonomic handle design makes flossing an easy, one handed operation.
- Comes with 15 flossing tips.
- Replaceable AA battery.



WaterPik®
Power
Flossers

treat
yourself
better



Includes :

- Flossing Tips (15)

WaterPik® WP-861 Complete Care 5.0

Why just brush when you can get Complete Care. Brushing alone is not enough, even with a high-end Sonic Toothbrush. WaterPik® Complete Care 5.0 combines Water Flosser and Sonic Toothbrush technologies to create a perfect daily oral care routine.



Brushing and flossing is now possible with one convenient device that saves counter space and power outlets. And it's perfect for anyone with braces or implants, and those with crowns, bridges, and veneers.

Includes :

- 2 - Classic Jet Tip
- 1 - Orthodontic Tip
- 1 - Plaque Seeker® Tip
- 1 - Pik Pocket™ Tip
- 2 - Triple Sonic/ Complete Care 5.0 Replacement Brush Heads - White

Reservoir capacity of 400 mL (15 oz)



Dentos India Pvt. Ltd.

301, Lavlesh Court, Pandit Varde Road, Bandra (W),
Mumbai - 400 050. • Tel.: 4250 5100 / 2642 1688
Fax : 4250 5101 • Email : mail@dentosindia.com
Website : www.dentosindia.com

NEW

The New Generation of

WaterPik®

Oral Irrigators
for Healthy,
Resilient Gums



gingivital
system